# TPH Spectra™

## Universal Composite Restorative

DENTSPLY

### DIRECTIONS FOR USE – ENGLISH

For dental use only.  
USA: Rx only.

#### 1. PRODUCT DESCRIPTION

TPH Spectra™ Universal Composite Restorative is a visible light cured, radiopaque, composite restorative material for anterior and posterior restorations and cosmetic veneering. It is to be used following the application of a suitable dentin/enamel adhesive and is compatible with all DENTSPLY adhesives designed for use with visible light cured composite restoratives (see complete Directions for Use of selected adhesive). Use of other dentin/enamel adhesive systems with TPH Spectra™ Composite Material is at the discretion and sole responsibility of the dental practitioner.

##### 1.1 Delivery forms\*

TPH Spectra™ Composite material is available in:

- Low Viscosity and High Viscosity
- Predosed Compules® Tips
- 3g Easy•Twist syringes

\*Some delivery forms may not be available in all countries.

Both Low and High Viscosity TPH Spectra™ Composite materials are available in the most commonly used Vita® shades. Additional complementary shades are available in Low Viscosity for esthetically demanding applications. Refer to the available shade conversion table for matching Prisma® and Bioform® shades to available equivalent to TPH Spectra™ Composite material shades.

TPH Spectra™ Composite material may be used with separately available methacrylate-based tints and opaques at the discretion and responsibility of the dental practitioner for individual custom characterization of the ultimate esthetic direct-placement restoration.

##### Basic Vita® Shades – Available in HV and LV

###### Light Shades (2)

Light Reddish Brown .....A1  
Light .....B1

###### Universal Body Shades: Reddish Brown (3)

Light Medium .....A2  
Medium .....A3  
Dark Medium .....A3.5

###### Universal Body Shades: Gray (1)

Medium .....C2

###### Characterization Shades (1)

Reddish Brown .....A4

##### Complementary Shades – Available in LV CUSTOM SHADES

Description	TPH Spectra™ Shade	VITA® SHADES Description	TPH Spectra™ Shade
<b>Light Shades (5)</b>			
Bleach White.....	BW		
Extra Light .....	XL		
Light.....	L		
Light Gray .....	LG	Light Gray.....	C1
<b>Universal Body Shades: Reddish Brown (2)</b>			
Light Medium.....	LB		
Dark Medium .....	DB		
<b>Universal Body Shades: Reddish Yellow (2)</b>			
		Light Medium .....	B2
		Dark Medium.....	B3
<b>Universal Body Shades: Gray (2)</b>			
Medium.....	LYG	Dark Medium.....	C3
<b>Characterization Shades (5)</b>			
Dark Reddish Brown.....	DY		
Medium Gray Opaque.....	C2-O	Gray .....	C4
Gray Brown.....	XGB	Reddish Gray.....	D3
<b>Enamels (3)</b>			
Clear Translucent Enamel.....	CE		
Yellow Translucent Enamel .....	YE		
Light Incisal.....	B1-I		

##### 1.2 Composition

Urethane modified Bis-GMA resin; TEGDMA Triethyleneglycoldimethacrylate; Polymerizable Dimethacrylate Resin; Camphorquinone (CQ) photoinitiator; Ethyl-4(dimethylamino)benzoate photoaccelerator; Butylated hydroxy toluene (BHT); UV stabilizer; Silanated barium-alumino-borosilicate glass; Silanated barium-boron-fluoro-alumino-silicate glass; Silicon dioxide; Fluorescent agent; Synthetic Inorganic Iron oxide pigments; Titanium dioxide

##### 1.3 Indications

1. TPH Spectra™ Composite material is indicated as a direct restorative for all cavity classes in anterior and posterior teeth.
2. TPH Spectra™ Composite material is indicated as an anterior direct esthetic veneering material.
3. TPH Spectra™ Composite material is indicated for the indirect fabrication of inlays and onlays.

##### 1.4 Contraindications

TPH Spectra™ Composite material is contraindicated for use with patients who have a known hypersensitivity to methacrylate resins.

##### 1.5 Compatible adhesives

TPH Spectra™ Composite material is used following application of a suitable dentin/enamel adhesive and is chemically compatible with DENTSPLY methacrylate-based dentin/enamel adhesives designed for use with visible light cured composite restoratives (see complete Directions for Use of selected adhesive).

### 2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions for use.

##### Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

##### 2.1 Warnings

TPH Spectra™ Composite material contains polymerizable methacrylates which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

##### 2.2 Precautions

1. This product is intended to be used only as specifically outlined in these Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and is the sole responsibility of the practitioner.
2. Contact with saliva and blood during composite placement may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
3. Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
4. Material should extrude easily DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause the Compules® Tip to rupture or to eject from the Compules® Tips Gun.
5. Compules® Tips may be used for direct intraoral application of restorative material into a cavity or for indirect application by first placing the restorative material on a pad. After intraoral use, discard Compules® Tips and do not reuse in other patients in order to prevent cross-contamination.
6. The TPH Spectra™ Composite material Easy•Twist syringe should be tightly closed immediately after use.
7. Use of Compules® Tips with the DENTSPLY Compules® Tips Gun is recommended.
8. Interactions:
  - Eugenol and hydrogen peroxide containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
  - TPH Spectra™ Composite material is light-cured material. Therefore, it should be protected from ambient light. Proceed immediately once the material has been placed.
  - If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.

### 2.3 Adverse reactions

Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage. Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous membranes: inflammation, edema, sloughing. See Warnings.

### 2.4 Storage

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Keep out of direct sunlight and store in a well ventilated place at temperatures between 2-24°C/35-75°F. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

## 3. STEP-BY-STEP INSTRUCTIONS

### Direct Restoration – Cavity Restoration, Cosmetic Veneering

#### 1. Shade Selection

Before selecting the shade, teeth should be clean, hydrated and free of extrinsic material or stain. The TPH Spectra™ Composite shade guide is to be used to accurately demonstrate the shade of the composite material. The TPH Spectra™ Composite shade guide tab is manufactured from the TPH Spectra™ Composite material of the selected shade.

TPH Spectra™ Composite has been scientifically designed to allow the tooth structure enamel and dentin to blend together with the composite with a life-like translucent result. The basic shades of TPH Spectra™ Composite adequately reproduce many of the traditional Vita® shades. See the shade coverage table below:

	B1	A1	A2	C2	A3	A3.5	A4
B1	●	●	●	●	●	●	●
A1	●	●	●	●	●	●	●
B2	●	●	●	●	●	●	●
D2			●	●	●	●	●
A2			●	●	●	●	●
C1			●	●	●	●	●
C2			●	●	●	●	●
D4			●	●	●	●	●
A3			●	●	●	●	●
D3			●	●	●	●	●
B3			●	●	●	●	●
A3.5			●	●	●	●	●
B4			●	●	●	●	●
C3			●	●	●	●	●
A4			●	●	●	●	●
C4			●	●	●	●	●

Additional Complementary shades are available for extended characterizations.

#### 2. Cavity Preparation

- Prepare the cavity so that no residual amalgam or restorative material is left.
- Rinse surface with water spray and carefully dry it with air spray. Do not desiccate the tooth structure.
- Use a dental dam or cotton rolls to isolate the cavity from contamination.

#### 3. Placement of Matrix

The use of a Mylar, sectional (e.g. Palodent®/Palodent® Plus Sectional Matrix System) or thin matrix band (e.g. AutoMatrix® Retainerless Matrix System) and subsequent burnishing of the matrix band will improve final interproximal contact and contour. Pre-wedging/ring placement is recommended.

#### 4. Pulp protection, tooth conditioning/dentin pretreatment, adhesive application

Refer to adhesive manufacturer's directions for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of TPH Spectra™ Composite material.

#### 5. Placement of TPH Spectra™ Universal Composite Restorative

##### Easy•Twist Syringe:

- Remove the cap.
- Turn the handle of the syringe slowly in a clockwise direction and dispense the necessary amount of the material onto a mixing pad.
- Point the front tip of the syringe upwards and turn the handle anti-clockwise to prevent oozing of the material.
- Immediately close the syringe with the cap.
- Protect the restorative material on the mixing pad against light.

##### Pre-dosed Compules® Tips:

#### Danger of injury due to excessive force

- Apply slow and steady pressure on the applicator gun
- Do not use excessive force – Compules® Tip rupture or ejection from applicator gun may result.

- Insert a Compules® Tip into the notched opening of the Compules® Tips Gun. Be certain that the collar on the Compules® Tip is inserted first.
- Remove the colored cap from the Compules® Tip. The Compules® Tip may be rotated 360° to gain the proper angle of entrance into the cavity or to the mixing pad. To dispense the material into a cavity preparation, use a slow, steady pressure. Do not use excessive force.
- To remove the used Compules® Tip, be sure that the Compules® Tips Gun plunger is pulled back completely by allowing the handle to open to its widest position. Apply a downward motion to the front end of the Compules® Tip and remove.

#### 5.1 Optional Flowable Liner

Placement of a compatible flowable liner (available separately) prior to placement of TPH Spectra™ Composite material is optional. Follow manufacturer's directions for use.

Dispense TPH Spectra™ Composite material directly into the cavity preparation/tooth surface from the Compules® Tip using slow, steady pressure. Excessive force is not necessary. Alternatively, material may be expressed onto a clean pad from the Compules® Tip or Easy•Twist syringe and carried to the preparation with suitable placement instrument. Adapt, contour and shape with appropriate composite instruments. Ensure material is intimately adapted to cavity walls prior to curing. Material may be placed and light cured in increments up to 2mm (see Curing, Step 6). **Technique Tip:** For esthetic blending of shades, shade layers may be individually cured to form a foundation for subsequent shade placement. TPH Spectra™ Composite material resists slumping, allowing carving of the majority of anatomical form prior to visible light curing.

#### 6. Curing

Light cure each area of the restoration surface with a suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470nm. Minimum light output must be at least 550mW/cm² exposure for at least 20 seconds at a distance not greater than 5mm from surface. Some advanced performance curing units have been shown to cure 2mm increments of most shades of TPH Spectra™ Composite material in 10 seconds. Refer to curing light manufacturer's recommendations for compatibility and curing recommendations. The TPH Spectra™ Composite material should be additionally exposed to the curing unit through the proximal, lingual, and buccal enamel walls following matrix removal for the recommended time.

#### Inadequate polymerization due to insufficient curing

- Check compatibility of curing light.
- Check curing cycle.
- Check curing output before each procedure.

#### 7. Finishing and Polishing

- Contour the restoration using finishing burs or diamond finishing instruments. Additional finishing is recommended by the use of Enhance® Finishing System. See manufacturer's complete directions for use.
- To achieve a very high luster on TPH Spectra™ Composite material, it is necessary to complete the polishing. PoGo® One Step Diamond Micro-Polisher System and/or Prisma®Gloss™ Composite Polishing Pastes are recommended. See manufacturer's complete directions for use.

### Indirect Restoration – Inlay, Onlay

#### 1. Shade Selection

Final desired shade selection should be accomplished prior to tooth preparation. Refer to the Shade Selection section above.

#### 2. Preparation

Design requirements are essentially a conventional preparation. Rounding of internal angles and refinement of the cavo-surface margin for enhancement of enamel bonding and finishing procedures is recommended. Refer to adhesive and/or luting cement manufacturer's directions for base/liner/pulp protection requirements.

#### 3. Impression, Master Cast Fabrication

Make an accurate impression of the preparation, prepare master cast and articulate opposing cast per usual technique

#### 4. Restoration Fabrication

- Block out any preparation undercuts present. Apply separating medium and/or die spacer to separated die. In most cases, the restoration may be fabricated in no more than 3 increments, each up to 2mm depth.
- Apply first increment/shade, creating restoration body. Adapt to die, staying short of margins. Recommended light curing is accomplished by placing in Triad® Light Curing Unit (DENTSPLY Trubyte) for 2 minutes. Apply second, body layer, allowing cut-back for occlusal anatomy and final proximal contour. Repeat Triad® Light Curing. Prior to final "enamel" shade placement, custom staining with appropriate compatible stains may be accomplished per manufacturer's instructions. Apply final layer, slightly overfilling and covering all margins. Replace die into articulated model. Establish all external, proximal and occlusal contacts and anatomy. Slight lubrication of adjacent and opposing model teeth is recommended. Cure briefly (10 seconds) with hand-held visible curing light unit to fix contours. Remove die with restoration, place in Triad Unit for final 2 minute curing.
- When removing restoration from die, it may be necessary to scrape die stone away from restoration margins, to prevent accidental chipping of restoration. Clean any residual die stone from restoration. Carefully trim away visible flash beyond preparation margins with acrylic bur.
- Gently sandblast the internal surfaces with 50m alumina abrasive. Remove visibly undercut material.
- Seat restoration onto preparation of uncut master cast, making adjustments as needed. Check for marginal integrity and overall fit and contour. Additional increments may be added if needed by roughening surface and applying compatible adhesive per manufacturer's instructions, followed by placement and light curing as outlined above.

#### 5. Finishing and Polishing (laboratory)

Complete as outlined in the above Finishing and Polishing Section for Direct Restorations.

#### 6. Cementation

Recommended cementation technique is bonding using a dual cure adhesive and esthetic resin cement. Follow adhesive and cement manufacturer's instructions for prepared tooth and restoration surface pretreatments.

#### 7. Adjusting, Finishing and Polishing (clinically)

Following cementation, make all necessary occlusal adjustments and polish any clinically adjusted surfaces as outlined in the above Finishing and Polishing Section for Direct Restorations.

#### 4. HYGIENE

##### Cross-contamination of Compules® Tips

- After intraoral use do not reuse Compules® Tips.
- Properly dispose the intraorally used and/or contaminated Compules® Tips in accordance with local regulations.

##### 4.1 Cleaning and disinfection of the syringe

To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid package contamination. Disinfect contaminated syringes with a water-based hospital-level disinfection solution according to national/local regulations.

Repeated disinfection may damage label. **NOTE:** Vigorous wiping can destroy the label. Wipe syringe gently.

##### 4.2 Cleaning and disinfection of TPH Spectra™ Shade Guide

- Disinfect the shade guide holder and individual tabs with a water-based hospital-level disinfection solution according to national/local regulations. Use of chlorine- (bleach), glutaraldehyde-, phenolic-, iodophor- or organic solvent (e.g. alcohol) based or containing products may over time affect the color shade.
- During disinfection, remove individual tabs from the shade guide holder.
- Follow the disinfectant manufacturer's Instructions for Use.
- After disinfection, clean the shade guide holder and individual tabs by scrubbing with water and soap.
- Rinse and dry.
- Do not autoclave the shade guide holder and individual tabs.

#### 4.3 Sterilizing the Compules® Tips Gun

##### Disassembling

- Partially close the applicator gun and place your thumb on the rear part of the hinge.
- Push upwards and lift hinge. The applicator gun is separated in two parts and the plunger is exposed.
- Remove restorative material with a soft paper tissue and alcohol of 70%.
- The Compules® Tips Gun may be cleaned by scrubbing with hot water and soap or detergent

##### Sterilization

Sterilize the applicator gun by steam autoclave (2.1-2.4 bar and 135-138°C, 30-35 psi and 275-280°F) 3 minutes minimum holding time. It is not recommended to submerge Compules® Tips Gun into disinfection solutions.

##### Reassembling

- After sterilization insert plunger into the applicator gun cylinder.
- Press components together and snap hinge mechanism in place.
- Do not continue to use damaged or soiled applicator guns.

#### 5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard uses: "YYYY/MM"
2. Compules® Tips Gun is warranted for one year from date of purchase.
3. The following numbers should be quoted in all correspondence:
  - Reorder Number
  - Lot number
  - Expiration date

Vita is not a registered trademark of DENTSPLY International.

Manufactured by:  
DENTSPLY Caulk  
38 West Clarke Avenue  
Milford, DE 19963 USA  
Tel.: 1-302-422-4511  
www.dentsply.com

DENTSPLY DE TREY GmbH  
De-Trey-Str. 1  
78467 Konstanz  
Germany  
Tel.: 49-7531-583-0  
www.dentsply.de

Distributed by:  
DENTSPLY Canada  
161 Vinyl Court  
Woodbridge, Ontario  
L4L 4A3 Canada

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DENTSPLY DeTrey GmbH